



Software Patent Institute Database of Software Technologies

Record Display

[New Search](#)[Modify Search](#)[Search Results](#)

Record 1

Serial number TDB1282.0113

Field Name	Contents of Record 1
Size of Record	2320 total bytes in record, 2108 in TX field
Title	Laser Activated Cleaning and Etching System and Method. December 1982.
Publication Date	December, 1982
Copyright Notices	<p>This record was retrieved 15:14:35 Sat, Jan 22, 19100 from the SPI Database of Software Technologies Copyright © 1995-1998 Software Patent Institute.</p> <p>Database entry Copyright © Software Patent Institute This article is © Copyright IBM Corp.</p>
Text of Submission	<p>2p. This is a system and method of treatment of a surface to be inspected or coated or otherwise treated in a vacuum system. This system provides a laser -activated gas which reacts with the surface of a sample to clean the surface or to thin the surface. This system can be employed to treat a surface located within a surface analyzer, a vacuum system or an electron microscope.</p> <p>The sample is placed within the apparatus which is to operate upon it. Located within the apparatus is a movable nozzle which supplies an inert gas which is capable of being energized by laser radiation to react locally with the surface of the sample. The gas reacts locally with the surface in the regions irradiated by the laser. For example, an inert polyatomic gas such as SF(6), reacts with Si in response to a focussed argon ion laser. The gases emitted during the reaction of the gas with Si are inert. Thus, etching of the Si surface can be performed rapidly without producing contamination. As soon as the inert gases have been evacuated from the system, the apparatus can be operated to produce the measurement or other result required.</p> <p>Fig. 1 shows a TEM (transmission electron microscope) system employing this system and method. A vacuum column extends vertically housing electron lenses B. A specimen holder C extends through the side wall of column A. A specimen D is supported by C. A movable nozzle E is aimed to direct a gas at the specimen D. The lower wall of the microscope includes a viewing port F through which a laser beam G is directed to reflect from a movable mirror H to strike specimen D on predetermined locations on its surface.</p>

Fig. 2 shows a vacuum deposition system employing this method. In this case **laser** beam J is directed through window K in the wall W of the vacuum chamber 0 to strike mirror M within the chamber. The **laser** beam Q is reflected from mirror M onto substrate L. A movable nozzle N is supported by a feedthrough R in opposing wall X of the vacuum chamber 0. The evaporation sources or sputtering targets are located in direction P below the substrate L.

Reference (pointer to work)	IBM TDB 12-82 p3775-3776 Order: 82A 62055
Submission Date	April, 1993
Date Loaded into Database	February 14, 1997
Publisher	IBM Corporation.
Journal	IBM TDB
Corporate Source	IBM
Country of Origin	U.S.A.
Publication Language	English
Source Type (Journal, book, etc.)	journal

[New Search](#)[Modify Search](#)[Search Results](#)[Help \(Manual pages\)](#)[Search the Manual](#)[SPI Home](#)[Comments](#)

This page was last revised Jan 14, 2000 by wem, bill@spi.org.